

Габаритные размеры																							
ТИП	DNA	DNM	a	g	f	x1	x2	n1	n2	n3	h1	h2	h3	m1	z1	s	b	i1	i	L	d	t	u
50-80/2	80	50	199	269,5	383	458	428	280	212	160	180	180	180	60	60	15	60	120	109	80	38	41,3	10
50-80/3	80	50	277	269,5	383	536	506	280	212	160	180	180	180	60	60	15	60	120	109	80	38	41,3	10
50-80/4	80	50	355	269,5	383	614	584	280	212	160	180	180	180	60	60	15	60	120	109	80	38	41,3	10
50-80/5	80	50	433	269,5	383	692	662	280	212	160	180	180	180	60	60	15	60	120	109	80	38	41,3	10
50-80/6	80	50	511	269,5	383	770	740	280	212	160	180	180	180	60	60	15	60	120	109	80	38	41,3	10
50-80/7	80	50	589	269,5	383	848	818	280	212	160	180	180	180	60	60	15	60	120	109	80	38	41,3	10
50-80/8	80	50	667	269,5	383	926	896	280	212	160	180	180	180	60	60	15	60	120	109	80	38	41,3	10
50-80/9	80	50	745	269,5	383	1004	974	280	212	160	180	180	180	60	60	15	60	120	109	80	38	41,3	10
50-80/10	80	50	823	269,5	383	1082	1052	280	212	160	180	180	180	60	60	15	60	120	109	80	38	41,3	10
50-80/11	80	50	901	269,5	383	1160	1130	280	212	160	180	180	180	60	60	15	60	120	109	80	38	41,3	10
50-80/12	80	50	979	269,5	383	1238	1208	280	212	160	180	180	180	60	60	15	60	120	109	80	38	41,3	10
50-80/13	80	50	1057	269,5	383	1316	1286	280	212	160	180	180	180	60	60	15	60	120	109	80	38	41,3	10

